

Drying Cabinets Features & Benefits

As you consider your options in selecting a drying cabinet that best suits your needs, the following points of difference between the ThermaSure dryers and others on the marketplace may help guide an informed decision.

◆ **ThermaSure dryers have the ability to hold on average 60% more circuits, hoses and tubes than competitive offerings.** In our 300 and 1000 Series models, the tube drying rack pulls out toward you making it easier and safer to load and unload these devices. The height of this rack is within easy reach, yet unlike the competitor's, it is high enough so that long circuits do not touch the floor.

◆ **The control panel is microprocessor controlled and is easy to use.** Our control panel provides you with the option of being able to select your own temperature, which will help to maximize the life of your medical devices. Your users merely program the exact drying cycle they want, or set the cycle to continuous mode, push a button and walk away. If you choose to configure your dryer with a HEPA filter, the control panel will automatically alert you when the filter needs to be changed.

◆ **The design of our ThermaSure dryers is more ergonomically correct than others on the market place.** There are two points of difference that set us apart and make us a better choice:

1. Our shelves are installed at heights that minimize bending and lifting for the user, which is especially important if you plan to dry heavy items or surgical trays.
2. Our filter system is located in an easy to reach area without the need to stand on a step stool or ladder and is accessible without the use of tools.

◆ **Robust stainless steel shelves provide maximum use of space and stand up to the test of time.** Our standard half-width shelves are tested to bear 40 pounds and optional full-width shelves will hold twice that. Only ThermaSure dryers can hold large trays, such as those that measure more than 18 inches in length. If you are no longer reprocessing hoses and tubes, you can easily configure the dryer with full-width shelves, maximizing flexibility and gaining broader application.

◆ **ThermaSure dryers are constructed of stainless steel for all parts that might come in contact with moisture (shelves, trays, floor, rivets), not aluminum.** While neither will rust, our cabinets will not dent or scratch as easily, and the heavier construction will hold up better over the system's lifespan. Aluminum is prone to corrosion.

◆ **Our current models of ThermaSure dryers may have as many as four fans distributed throughout the cabinet, depending on options chosen.** The placement of these heaters and fans is important in creating greater cfm (cubic foot per minute) airflow. Our drying cabinets are designed to evenly distribute air flow and heat throughout the cabinet, minimizing drying time and maximizing efficiency.

Premium Medical Device Drying Cabinets