


Model 540 Washer/Pasteurizer
Technical Specifications



PHYSICAL SPECIFICATIONS	
WIDTH:	26 1/2" (68 cm)
DEPTH:	31" (79 cm)
HEIGHT:	44" (112 cm)
CLEARANCE:	
(LID RAISED)	61" (155 cm)
WEIGHT (EMPTY):	205 lbs (93 kg)
(FULL):	413 lbs (187 kg)
VOLUME:	25 Gal. (95 L)

ELECTRIC SPECIFICATIONS	
VOLTAGE:	115 VOLTS (230 VOLTS available)
CYCLES:	60 HERTZ
CURRENT:	20 AMPS
REQUIRED CIRCUITS:	
	1) 120V/30A, LOCKING RECEPTACLE. (NEMA L5-30R)
	2) HIGH THROUGHPUT HEATER: SEE SPECS
	*EQUIPPED WITH 5' CORD AND 30 AMP LOCKING PLUG (NEMA L5-30P)

PLUMBING SPECIFICATIONS
WATER REQUIREMENTS:
3/4" HOSE BIB FROM HOSPITAL SUPPLY (120°F RECOMMENDED)
165°F FROM HIGH THROUGHPUT HEATER
DRAIN REQUIREMENTS: 1 1/2" ID (3.8 cm)

#224861 

High Level Disinfection Systems for Healthcare

Cenorin, LLC • 6324 South 199th Place • Suite 107 • Kent, WA 98032
 Phone: 253-395-2400 • Fax: 253-395-2650 • Toll Free: 1-800-426-1042