


**Model 540 230V Washer/Pasteurizer
Technical Specifications**



PHYSICAL SPECIFICATIONS		
WIDTH:	26 1/2"	(68 cm)
DEPTH:	31"	(79 cm)
HEIGHT:	44"	(112 cm)
CLEARANCE: (LID RAISED)	61"	(155 cm)
WEIGHT (EMPTY):	205 lbs	(93 kg)
(FULL):	413 lbs	(187 kg)
VOLUME:	25 Gal.	(95 L)

ELECTRIC SPECIFICATIONS	
VOLTAGE:	230 VOLTS
CYCLES:	50/60 HERTZ
CURRENT:	12 AMPS
REQUIRED CIRCUITS:	
	1) 240V/30A, LOCKING RECEPTICLE. (NEMA L5-30R)
	2) HIGH THROUGHPUT HEATER: SEE SPECS
*EQUIPPED WITH 5' CORD AND 30 AMP LOCKING PLUG (NEMA L5-30P)	

PLUMBING SPECIFICATIONS
WATER REQUIREMENTS:
3/4" HOSE BIB FROM HOSPITAL SUPPLY (120°F RECOMMENDED)
165°F FROM HIGH THROUGHPUT HEATER
DRAIN REQUIREMENTS: 1 1/2" ID (3.8 cm)

[High Level Disinfection Systems for Healthcare](#)

Cenorin, LLC • 6324 South 199th Place • Suite 107 • Kent, WA 98032
Phone: 253-395-2400 • Fax: 253-395-2650 • Toll Free: 1-800-426-1042